

OMG!!!

Helen Thom BMSEDG#9

Caring Expert Quality



Objectives

- To try and understand anti-G
- How to distinguish between anti-D, C and G
- Why its important
- Some example scenarios
- Anonymous poll at the end
- Keep an eye out for the clues





So what is the G antigen?

- Part of the Rh Blood group System
- Rh12
- Clinically significant
- Enhanced by papain (enzyme treated cells)
- Present on red cells which carry the D antigen or the C antigen or both









https://www.bbguy.org/2016/06/17/want-g-wiz/



Some common G pos and G neg phenotypes

Name	Pheno	G Pos	Name	Pheno	G neg
R1R1	DCe/DCe	D+ C+ G+	rr	dce/dce	D- C- G-
R2R2	DcE/DcE	D+ C- G+	r"r"	dcE/dcE	D- C- G-
r'r	dCe/dCe	D- C+ G+	r"r	dcE/dce	D- C- G-
R1R2	DCe/DcE	D+ C+ G+			
RO	Dce/Dce	D+ C- G+		Not reall	y 'common'
R1r	DCe/dce	D+ C+ G+			
R2r	DcE/dce	D+ C- G+			





Cell	Rh	D	С	E	с	е	м	N	S	S	P1	Luª	к	k	Kpª	Le ^a	Le ^b	Fyª	Fy⁵	Jkª	Jk⁵	Other	ΙΑΤ	EIAT
1	R ₁ ^w R ₁	+	+	0	0	+	+	0	+	+	0	0	0	+	0	0	+	+	0	+	0		3	5
2	R ₁ R ₁	+	+	0	0	+	0	+	0	+	0	0	+	+	0	+	0	0	+	0	+		3	5
3	R ₂ R ₂	+	0	+	+	0	+	0	+	0	0	0	0	+	0	+	0	0	+	+	0		4	5
4	r'r	0	+	0	+	+	0	+	+	0	0	0	0	+	0	0	+	0	+	+	0		0	0
5	r"r	0	0	+	+	+	0	+	0	+	4	0	0	+	0	0	+	+	0	0	+		0	0
6	rr	0	0	0	+	+	+	0	+	0	4	0	+	0	0	0	+	+	0	0	+		0	0
7	rr	0	0	0	+	+	0	+	0	+	2	+	+	+	0	+	0	0	+	+	0		0	0
8	rr	0	0	0	+	+	0	+	0	+	0	0	0	+	+	0	+	0	+	0	+		0	0
9	rr	0	0	0	+	+	+	0	0	+	2	0	0	+	0	+	0	+	0	0	+		0	0
10	rr	0	0	0	+	+	0	+	0	+	3	0	0	+	0	+	0	+	0	+	0		0	0
																						Auto	0	/
																						K control	2	/





Cell	Rh	D	с	E	с	е	М	N	S	S	P1	Luª	к	k	Kpª	Leª	Le ^b	Fyª	Fy⁵	Jkª	JkÞ	Other	ΙΑΤ	EIAT
1	$R_1^w R_1$	+	+	0	0	+	+	0	+	+	0	0	0	+	0	0	+	+	0	+	0		3	5
2	R ₁ R ₁	+	+	0	0	+	0	+	0	+	0	0	+	+	0	+	0	0	+	0	+		3	5
3	R_2R_2	+	0	+	+	0	+	0	+	0	0	0	0	+	0	+	0	0	+	+	0		0	0
4	r'r	0	+	0	+	+	0	+	+	0	0	0	0	+	0	0	+	0	+	+	0		3	5
5	r"r	0	0	+	+	+	0	+	0	+	4	0	0	+	0	0	+	+	0	0	+		0	0
6	rr	0	0	0	+	+	+	0	+	0	4	0	+	0	0	0	+	+	0	0	+		0	0
7	rr	0	0	0	+	+	0	+	0	+	2	+	+	+	0	+	0	0	+	+	0		0	0
8	rr	0	0	0	+	+	0	+	0	+	0	0	0	+	+	0	+	0	+	0	+		0	0
9	rr	0	0	0	+	+	+	0	0	+	2	0	0	+	0	+	0	+	0	0	+		0	0
10	rr	0	0	0	+	+	0	+	0	+	3	0	0	+	0	+	0	+	0	+	0		0	0
																						Auto	0	/
																						K control	2	/





Cell	Rh	D	С	E	с	е	м	N	S	S	P1	Luª	к	k	Kpª	Leª	Le ^b	Fyª	Fy⊳	Jkª	Jk⁵	Other	ΙΑΤ	EIAT
1	R ₁ ^w R ₁	+	+	0	0	+	+	0	+	+	0	0	0	+	0	0	+	+	0	+	0		3	5
2	R ₁ R ₁	+	+	0	0	+	0	+	0	+	0	0	+	+	0	+	0	0	+	0	+		3	5
3	R ₂ R ₂	+	0	+	+	0	+	0	+	0	0	0	0	+	0	+	0	0	+	+	0		4	5
4	r'r	0	+	0	+	+	0	+	+	0	0	0	0	+	0	0	+	0	+	+	0		3	5
5	r"r	0	0	+	+	+	0	+	0	+	4	0	0	+	0	0	+	+	0	0	+		0	0
6	rr	0	0	0	+	+	+	0	+	0	4	0	+	0	0	0	+	+	0	0	+		0	0
7	rr	0	0	0	+	+	0	+	0	+	2	+	+	+	0	+	0	0	+	+	0		0	0
8	rr	0	0	0	+	+	0	+	0	+	0	0	0	+	+	0	+	0	+	0	+		0	0
9	rr	0	0	0	+	+	+	0	0	+	2	0	0	+	0	+	0	+	0	0	+		0	0
10	rr	0	0	0	+	+	0	+	0	+	3	0	0	+	0	+	0	+	0	+	0		0	0
																						Auto	0	/
																						K control	2	/

Anti-G (cannot excluded anti-D or C)

		NHS
Blood	and	Transplant

Cell	Rh	D	С	E	с	е	м	N	S	S	P1	Luª	к	k	Kpª	Leª	Le ^b	Fyª	Fy ^b	Jk ^a	Jk⊳	Other	ΙΑΤ	EIAT
1	$R_1^w R_1$	+	+	0	0	+	+	0	+	+	0	0	0	+	0	0	+	+	0	+	0		4	5
2	R ₁ R ₁	+	+	0	0	+	0	+	0	+	0	0	+	+	0	+	0	0	+	0	+		4	5
3	R ₂ R ₂	+	0	+	+	0	+	0	+	0	0	0	0	+	0	+	0	0	+	+	0		3	4
4	r'r	0	+	0	+	+	0	+	+	0	0	0	0	+	0	0	+	0	+	+	0		3	4
5	r''r	0	0	+	+	+	0	+	0	+	4	0	0	+	0	0	+	+	0	0	+		0	0
6	rr	0	0	0	+	+	+	0	+	0	4	0	+	0	0	0	+	+	0	0	+		0	0
7	rr	0	0	0	+	+	0	+	0	+	2	+	+	+	0	+	0	0	+	+	0		0	0
8	rr	0	0	0	+	+	0	+	0	+	0	0	0	+	+	0	+	0	+	0	+		0	0
9	rr	0	0	0	+	+	+	0	0	+	2	0	0	+	0	+	0	+	0	0	+		0	0
10	rr	0	0	0	+	+	0	+	0	+	3	0	0	+	0	+	0	+	0	+	0		0	0
																						Auto	0	/
																						K control	2	/

What makes us think its anti-G: Spot the difference!





R2R2 cells have more D antigen sites than R1R1 (steric hinderance) so with anti-D + C, you expect to see stronger reaction vs R2R2

Cell	Rh	D	С	E	с	е	м	N	S	S	P1	Luª	к	k	Kpª	Leª	Le ^b	Fyª	Fy⁵		IAT	EIAT
1	R ₁ ^w R ₁	+	+	0	0	+	+	0	+	+	0	0	0	+	0	0	+	+	0		4	5
2	R ₁ R ₁	+	+	0	0	+	0	+	0	+	0	0	+	+	0	+	0	0	+	H ₁ H ₁	4	5
3	R ₂ R ₂	+	0	+	+	0	+	0	+	0	0	0	0	+	0	+	0	0	+	D+C+ = G+	3	4
4	r'r	0	+	0	+	+	0	+	+	0	0	0	0	+	0	0	+	0	+	+ 0	3	4

Weaker reaction vs R2R2 cell *indicates* anti-G. It *does not* exclude anti-D and/or C

BSH Guidelines: Selection of issue and Information of issue and Transplant Issue of Red cells

- 7.8.3.D negative red cells should always be selected for:
 - D negative women of childbearing potential (<51 years)
 - D negative patients <18 years old
 - Patients who have formed immune anti-D, even if not currently detectable
 - Transfusion-dependant D negative adults
- 7.10.2. Antigen negative red cells should also be selected when a clinically significant antibody has previously been identified, but cannot be detected or identified in the current sample

 7.10.3.Patients with anti-D who are rr (ccddee) should receive rr (D- C- E-), K negative blood

 7.10.4.Patients with other Rh antibodies should be additionally matched for C, c, E and e in order to prevent further Rh alloimmunisation, provided this does not impede delivery of effective transfusion support

https://b-s-h.org.uk/guidelines/guidelines/pre-transfusion-compatibility-procedures-in-blood-transfusion-laboratories/



Why does G matter?

- For transfusion purposes in an adult above child bearing potential it doesn't
 - Select rr (D- C- G-)
- For children, and those of childbearing potential it does
- The key thing with anti-G is to exclude anti-D (and/or anti-C)
 - Blood selection
 - Prevent sensitisation
 - Eligibility for RAADP
 - Assess risk for HDFN
 - -For anti-D/G: need to quant neat and adsorbed plasma
 - -Can't test anti-C/G in isolation titre vs r'r

How to we do it?



NHS Blood and Transplant

Sequential adsorptions with papainised cells

Mix patients plasma with r'r – D- C+ G+

Mix patients plasma with R2R2/Ror - D+ C- G+

Incubate for 10 mins at 37°C Centrifuge Remove plasma – save each aliquot Repeat adsorptions using each of the absorbed plasma – r'r, R2R2/Ror

Antigen positive adsorptions cells remove the corresponding antibody from the patient's plasma

Think antibody/antigen interactions

What you end up with....







If you have anti-D, C and G



Cell	Rh	D	с	E	С	е	М	N	S	S	P1	Luª	к	k	Kpª	Leª	Le ^b	Fyª	Fy⁵	Jkª	JkÞ	Neat IAT	r'r IAT	Ror IAT
1	R ₁ ^w R ₁	+	+	0	0	+	+	0	+	+	0	0	0	+		r'	r pla	sma:			0	4	2	3
2	R ₁ R ₁	+	+	0	0	+	0	+	0	+	0	0	+	+		r an	remo	ved	2		+	4	2	3
3	R ₂ R ₂	+	0	+	+	0	+	0	+	0	0	0	0	+	lea	aves	anti	-D b	ehino	d 🗌	o	3	2	0
4	r'r	0	+	0	+	+	0	+	+	0	0	0	0	+	0	0	+	0	+	+	0	3	0	3
5	r''r	0	0	+	+	+	0	+	0	+	4	0	0	+		Ro	l or pla	l asma			+	0	0	0
6	rr	0	0	0	+	+	+	0	+	0	4	0	+	0		r	emo	ved			+	0	0	0
7	rr	0	0	0	+	+	0	+	0	+	2	+	+	+	lea	an aves	ti-D a anti	and (-C b	G ehina	d	0	0	0	0
8	rr	0	0	0	+	+	0	+	0	+	0	0	0	+	+	0	+	0	+	0	+	0	0	0
9	rr	0	0	0	+	+	+	0	0	+	2	0	0	+		Ne	at p	lasm	a		+	0	0	0
10	rr	0	0	0	+	+	0	+	0	+	3	0	0	+	Т	here		cont	ains	Γ	0	0	0	0
																ar			כ		Auto	0	/	/
																					K control	2	/	/





Cell	Rh	D	с	E	С	е	м	N	S	S	P1	Luª	к	k	Kpª	Le ^a	Le ^b	Fyª	Fy ^b	Jk ^a	Jk⋼	Neat IAT	r'r IAT	Ror IAT
1	R ₁ ^w R ₁	+	+	0	0	+	+	0	+	+	0	0	0	+		r'	r pla	sma:			0	4	4	0
2	R ₁ R ₁	+	+	0	0	+	0	+	0	+	0	0	+	+		r an	emo ti-C a	ved	G		÷	4	4	0
3	R_2R_2	+	0	+	+	0	+	0	+	0	0	0	0	+	lea	aves	anti	-D b	ehino	d	o	3	3	0
4	r'r	0	+	0	+	+	0	+	+	0	0	0	0	+	0	0	+	0	+	+	0	3	0	0
5	r''r	0	0	+	+	+	0	+	0	+	4	0	0	+		Ro	br pla	asma	1:		+	0	0	0
6	rr	0	0	0	+	+	+	0	+	0	4	0	+	0		r		ved			+	0	0	0
7	rr	0	0	0	+	+	0	+	0	+	2	+	+	+	N	an Iothii	ng le	ft be	ے hind		0	0	0	0
8	rr	0	0	0	+	+	0	+	0	+	0	0	0	+	Т	here	fore	does	s not		+	0	0	0
9	rr	0	0	0	+	+	+	0	0	+	2	0	0	+	v	cor					+	0	0	0
10	rr	0	0	0	+	+	0	+	0	+	3	0	0	+		Ne	at pl	asm	а		0	0	0	0
															Т	here	fore		ains		Auto	0	/	/
																c					K control	2	/	/



If you have anti-C and G



Cell	Rh	D	С	E	с	е	М	N	S	s	P1	Lu ^a	к	k	Kpª	Le ^a	Le ^b	Fyª	Fy⁵	Jkª	Jk♭	Neat IAT	r'r IAT	Ror IAT
1	R ₁ ^w R ₁	+	+	0	0	+	+	0	+	+	0	0	0	+		r'ı	r pla:	sma:			0	4	0	3
2	R ₁ R ₁	+	+	0	0	+	0	+	0	+	0	0	+	+		r ant		ved	2		+	4	0	3
3	R ₂ R ₂	+	0	+	+	0	+	0	+	0	0	0	0	+	N	othir	ng le	ft be	5 hind		0	3	0	0
4	r'r	0	+	0	+	+	0	+	+	0	0	0	0	+	Т	here	fore	does	s not	:	0	3	0	3
5	r"r	0	0	+	+	+	0	+	0	+	4	0	0	+	v	cor	itain	anti-	·D		+	0	0	0
6	rr	0	0	0	+	+	+	0	+	0	4	0	+	0		Rc	or pla	asma	a:		+	0	0	0
7	rr	0	0	0	+	+	0	+	0	+	2	+	+	+		r an	emo ti-D a	ved	3		0	0	0	0
8	rr	0	0	0	+	+	0	+	0	+	0	0	0	+	Le	aves	beh	ind a	anti-	с [+	0	0	0
9	rr	0	0	0	+	+	+	0	0	+	2	0	0	+	0	No	ot pl				+	0	0	0
10	rr	0	0	0	+	+	0	+	0	+	3	0	0	+	т	here	fore	cont	a tains		0	0	0	0
																а	nti-C	C+G			Auto	0	/	/
																					K control	2	/	/





Cell	Rh	D	С	E	с	е	М	N	S	S	P1	Luª	к	k	Kpª	Leª	Le ^b	Fyª	Fy⁵	Jkª	Jkp	Neat IAT	r'r IAT	Ror IAT
1	R ₁ ^w R ₁	+	+	0	0	+	+	0	+	+	0	0	0	+		r'i r	r pla:	sma: ved			0	4	0	0
2	R_1R_1	+	+	0	0	+	0	+	0	+	0	0	+	+		an	ti- <mark>C</mark> a	and (G		+	4	0	о
3	R ₂ R ₂	+	0	+	+	0	+	0	+	0	0	0	0	+	N T	lothir	ng le	ft be	hind		O	3	0	0
4	r'r	0	+	0	+	+	0	+	+	0	0	0	0	+		cor	itain	anti-	·D		0	3	0	0
5	r"r	0	0	+	+	+	0	+	0	+	4	0	0	+	•				•		+	0	0	0
6	rr	0	0	0	+	+	+	0	+	0	4	0	+	0		RC r	or pia emo	asma ved	1:		+	0	0	0
7	rr	0	0	0	+	+	0	+	0	+	2	+	+	+		an	ti-D a	and (G		0	0	0	0
8	rr	0	0	0	+	+	0	+	0	+	0	0	0	+	N T	lothii here	ng le fore	ft be does	hind s not		+	0	0	0
9	rr	0	0	0	+	+	+	0	0	+	2	0	0	+		cor	tain	anti-	С		+	0	0	0
10	rr	0	0	0	+	+	0	+	0	+	3	o	0	+		Ne	at pl	asm	a		0	0	0	0
															Т	here	fore	cont	ains		Auto	0	/	/
																	anti	-G			K control	2	/	/



TYPICAL RESULTS

Absorbed with	Tes	st Cells	Conclusion
	r'	R2R2 / R _o	
r'	0	0	Anti-G only
R2R2 / R₀	0	0	
r'	0	+	Anti-D + G
R2R2 / R _o	0	0	
r'	0	0	Anti-C + G
R2R2 / R₀	+	0	
r'	0	+	Anti-D + C +/- G
R2R2 / R _o	+	0	



Case Scenario 1

- 65 year old male
- Hb 75
- 2 units RBC requested
- ABID: presents as apparent anti-D+C
- What would you do next?
 - Confirm anti-C using second cell
 - Select ABO and K compatible, D- C- E-





Case Scenario 2

- 7 year old female
- Multi transfused since birth
- No previous alloantibodies
- ABID: apparent anti-D+C
- Extended Rh phenotype mixed field vs anti-C well
 Unit was C+ (therefore G+)
- What would you do next?

YA	30	d.	A	15	018	41	42	43
			1	1	1			
	-	U	H					
Anti-C	Anti-c	Anti-E	Anti-e	Anti-K	Ctl	Anti-A	Anti-B	Anti-D VI
DP	++++		++++	•	-		-	

- G invest showed anti-C and G (n/s)
- Select ABO compatible, D- C- E- K-



Case Scenario 3

- Pregnant woman at 28 weeks' gestation
- RhD negative
- Second pregnancy
- Negative screen at booking
- ABID: apparent anti-D+C
- No RAADP
- Never transfused
- What would you do next? Refer to RCI



Case Scenario 3 continued

- Potential outcomes:
 - Anti-D + C quant anti-D and titre anti-C
 - Anti-D + C + G quant anti-D vs neat and adsorbed, titre anti-C/G
 - Anti-D + G quant anti-D vs neat and adsorbed, titre anti-G
 - Anti-C + G titre anti-C/G quant not required
 - Anti-G titre anti-G quant not required
- Quantitate neat and adsorbed plasma: the presence of anti-G can elevate quant result (R1R1 cells used) – need to assess risk of HDFN (using BSH guidelines) and eligibility for RAADP
- Can't titre C/G in isolation titre vs r'r and report combined result



Case Scenario 3 continued

Antibody Quantification

Anti-D concentration	Predicted clinical outcome
Less than 4 IU mL ⁻¹	HDFN unlikely, continue to monitor
4-15 IU mL ⁻¹	Moderate risk of HDFN, requiring referral to a fetal medicine specialist
More than 15 IU mL ⁻¹	High risk of HDFN requiring referral, as above

Antibody Titration
Anti-C
Anti-G
Anti-C/G – can't titre in isolation
r'r cell: D- C+ G+
 <32: low risk HDFN
 ≥32: high risk HDFN

https://b-s-h.org.uk/guidelines/guidelines/blood-grouping-and-antibody-testing-in-pregnancy/

Blood and Transplant Example of Quant Results for anti-D with G (+/-C)



- Neat plasma: 5.0 IU/mL
 - elevated due to the presence of anti-G
- r'r adsorbed plasma: <0.2 IU/mL
 - Removed anti-C and G leaves anti-D
 - N/S anti-D (likely prophylaxis)
 - Check history
- This patient is eligible for RAADP

- Neat plasma: 5.0 IU/mL
 - elevated due to the presence of anti-G
- r'r adsorbed plasma: 4.0IU/mL
 - Removed anti-C and G leaves anti-D
 - In this case, there is an alloanti-D
 - Moderate risk HDFN
- This patient is NOT eligible for RAADP

Propholicious?

Antibodies	Anti-D specifi using adsor	city and level bed plasma	Eligible for RAADP
Anti-D + C	N/S	<0.4 IU/mL	Yes
Anti-D + C	ALLO	Any level	No
Anti-D + C + G	N/S	0.4 IU/mL	Yes
Anti-D + C + G	ALLO	Any level	No
Anti-D + G	N/S	<0.4 IU/mL	Yes
Anti-D + G	ALLO	Any level	No
Anti-C + G	N/A	N/A	Yes
Anti-G	N/A	N/A	Yes

It isn't as clear cut as this, but it gives you an idea of the eligibility for RAADP The values obtained using the adsorbed plasma help to inform these decisions RCI medics and Clinical Scientists are available to make the clinical decisions



Recap of the Clues



Weaker reaction vs R2R2 cell indicates anti-G



Patient's phenotype: must be G negative (rr, r"r" or r"r)



Elevated quant result



You will probably have a headache!



Thank you for listening



• I hope this has helped you to understand anti-G

- Please take our poll
- All answers anonymised