



Referrals to RCI

Eve Hayes- Biomedical Scientist Red Cell Immunohaematology (RCI) – Colindale





Requesting tests




Requesting tests









Requesting Tests



- Be specific!
- If requesting a Weak D investigation or ABO subgrouping then state this!

 If you're not sure, then consider putting your results and state where the problem is!

Sending samples



FAO: RED CELL IMMUNOH	AEMATOLOGY
-----------------------	-------------------

NHS Blood and Transplant - Colindale

Charcot Road, Colindale London

ROUTINE / URGENT

NW9 5BG

Diagnostic Specimen

STORE BETWEEN 2°C - 8°C

- Clear labelling
- Specifies which department

- Urgent or Routine?
- NHSBT centre with RCI or Interim site?



https://hospital.blood.co.uk/diagnostic-services/red-cell-immunohaematology/labels/

1A Request Form

Blood and Transplant

Request for Reference Serology 1

Please see reverse for samples required Place correctly labelled samples in bag, remove protective strip, fold onto bag and seal firmly

At least three points of matching identification must be used on forms and sample tubes (see over)

Please do not use this form for routine antenatal requests Please write clearly inside text boxes, one character per box and mark option boxes as appropriate TEAR HERE PATIENT DETAILS First name DOR number **HOSPITAL & REQUESTER DETAILS** Phone Paternal or cord samples: Please give details of linked antenatal patient: CROSSMATCH REQUEST INVESTIGATIONS REQUESTED Date & time required required Special requirements: (e.g. CMV neg, irradiated) REFERRING LABORATORY'S FINDINGS YOUR RESULTS (state technique(s) & date tested) Antibody screen: Blood group: Ab ID: Sample date & time Other results: **DIAGNOSIS & CLINICAL HISTORY** Diagnosis: Hb & date: Previous transfusion? Yes Anti-D Ig (Dose & date given) Signature

NHSBT USE

Type of sample if not EDTA:

NHSBT sample number

Date / time sample received

Hospitals and Science website:

https://hospital.blood.co.uk/diagnostic-services /red-cell-immunohaematology/

Urgent Requests

RCI Record of Verbal Contacts - Information Received

Location				Phone number		
Hospital contact name (record full name	e)			Job title / Position		
PATIENT DETAILS						
Full Name						
Date of birth			Hospital N NHS num			
CLINICAL CONDITION	ON/TYPE OF PA	TIENT/SERO	LOGICAL F	INDINGS:		
Type of patient (ple	ase circle):	Medical	/ haematolo	gy / obstet	ric / surgical	
Date of last transfus	sion:					
Hb (Result, Time, Date):						
Blood group / Antib DAT result:	odies /					
Special requirement	ts (please circle)): CMV - / I	rradiated / V	Vashed / HI	b\$ - / Other / I	None
Additional informati (anti-D given, Consultant details, findings etc)	·					
SAMPLE TRANSPO	RT DETAILS:					
Transport method;		arrival o sample:				
TRANSFUSION REG	QUIREMENTS (p			ot an acce	ptable time/da	ite)
No. of units requested:		Date and required	d time I at hospital	:		
Transport of units (please circle): R	outine / Ad I	noc / Hospit	al collect /	Blue light / No	ot stated
Call Taken by (sign):		Date:			Time:	VMB updated? (circle) Y N/A



- Every time an urgent sample is called to RCI we have to complete one of these forms and an alternative form for FMH samples.
- Why?

RCI Techniques



- LISS
- NISS
- IAT
- DAT
- Neutralisation (Lewis, P1, Chido/Rogers, Knops)
- Adsorptions
- Elution
- Flow cytometry
- Quantification using continuous flow analyser (CFA)
- PCR

Abbreviations



- BRIAT Bio-Rad IAT
- EIAT- Bio-Rad IAT
- LISS or LIAT LISS IAT

- NAAD or NAD- No Alloantibodies Detected
- SNDT- Specificity Not Determined



Common Cases-1

- Patient not seen at this hospital before- no other clinical information!
- 51 year old patient
- Grouped as B RhD Positive, previously typed as R1r K-
- Automated panel on IH1000

Common Cases-1




Sample Name PATIENT					Reques	tor			_		Data	base R	ef No.		-		Tested by		*I	H1000			
Date of Birth			_			Hosp n	0			_		Sample No.				_		Date Tested Product			_		
			Lot No.		Produc	t			Lot No.		Produ	ct			Lot	No.	Lot No.						
Panel in Alsevers				R144 34	66	Panel in		b		R143 34	166	Panel in CellMedia					3 3466	Panel in LISP			R146 3466		
Panel P	apainised	d in Alsevers R154 3466			66		apainise		IIStab	R153 3466		Panel	Papani	sed in C	ellMedia	R17	3 3466	EXPIRY DATE: 2020.11.19					
	Rh	M	N	s	s	P1	Lua	К	k	Kpa	Lea	Le ^b	Fy ^a	Fy ^b	Jka	Jkb	Other	* BRIA	FIAT				
1	R ₁ wR ₁	0	+	0	+	3	0	0	+	0	0	+	0	+	0	+		1	1				
2	R ₁ R ₁	+	0	+	0	0	0	+	+	0	0	0	+	0	+	0		2	2				
3	R ₂ R ₂	+	0	0	+	2	0	0	+	0	+	0	+	0	0	+	Cob+, Bgb		2				
4	r'r	0	+	0	+	2	0	0	+	0	+	0	+	+	+	0		1	2				
5	r"r	0	+	+	0	1	0	0	+	0	0	+	0	+	+	0		2	2				
6	rr	+	+	0	+	1	0	+	0	0	0	+	0	+	+	0		1	2				
7	rr	0	+	0	+	1	0	+	+	0	+	0	0	+	0	+		2	2				
8	rr	0	+	0	+	2	0	0	+	+	0	+	+	0	+	0		1					
9	rr	+	0	+	0	0	0	0	+	0	+	0	+	0	0	+	Cob+	2					
10	rr	+	0	0	+	4	+	0	+	0	0	+	0	+	0	+		1					
Auto																		0					
											700 -37.50						7						
Antibo	dy Titre						•							Grou	p.			Pheno)				
ilution:												Cell	ld										
nti-														Conc	lusion								
rchive																							
nti																							
rchive																							
AT B	atch No							•															
PS		IgG		IgA		IgM		C3c		C30	1	Ct	I										
														Repor	rting BN	1S:		(Checkir	ng BMS			
Reagent	t			T						T			Γ					T		T			
atch										+													



Common Case 1

- Patient record from sample sent to RCI Tooting stated the patient was on Daratumumab.
- From these results we would suspect Daratumumab and check with the hospital if they were on the treatment due to results.



What testing should we do before starting treatment?





Common Case 2

- New patient
- 34 year old male
- Pre-op sample
- Anti-M stated on request form.
- Grouped by IH1000 as O RhD negative
- Rh/K- C-c+E-e+K-

Common Cases-2

Pipette batch No's



NHSBT Reagents Panel 1 **Antibody Investigation Worksheet** Sample Name PATIENT Tested by Requestor Database Ref No. *IH1000 Date of Birth **Date Tested** Hosp no Sample No. Product Lot No. Product Lot No. Product Lot No. Product Lot No. Panel in Alsevers R144 3466 Panel in CellStab R143 3466 Panel in CellMedia R163 3466 Panel in LISP R146 3466 R154 3466 Panel Papainised in Alsevers Panel Papainised in CellStab R153 3466 Panel Papanised in CellMedia R173 3466 EXPIRY DATE: 2020.11.19 BRIAT EIRT Rh N S Lua Kpa Lea Jka Jkb Other R₁^wR₁ 1 0 0 3 0 0 0 0 0 0 0 0 + + + + + 0 2 R₁R₁ 0 0 + 0 0 0 + 0 0 + 0 + 0 3 R_2R_2 0 0 2 0 0 0 0 0 Cob+, Bgb+ WF + + + + 0 + 4 r'r 0 0 2 0 + 0 0 0 0 0 5 r"r 0 0 + 0 1 0 0 + 0 0 0 0 + + 6 rr 0 1 0 0 0 0 0 0 0 + + 7 rr 0 + 0 1 0 0 0 0 0 0 + + + 8 0 0 0 0 + 0 + 0 + + 0 + 0 9 rr 0 0 0 0 0 0 0 0 0 Cob+ 10 rr 0 0 4 0 0 + + 0 0 0 + 0 Auto Group Pheno **Antibody Titre** Dilution Cell Id Conclusion Anti-Archive Anti Archive **DAT Batch No** PS **IgG** IgA lgM C3c C3d Ctl Reporting BMS: Checking BMS..... Reagent Batch No's

G:\009 DD&R\003 Red Cell Immunohaematology\001 Everyone\Reagent product profiles



Interesting Case 1

- 28 year old pregnant woman
- Two previous pregnancies
- First referred to RCI in 2015
- Grouped as B RhD Positive by IH1000
- Srilankan put as ethnicity on request form
- Previously detected a high frequency antibody
- EDD: May 2021

Interesting Case 1- Serology Blood and Transplant



	Reagen		anel							Antiboo	iy inve	stigati	on wo	orksne	et							
Date of Birth PATIENT I					Reques	tor		_	-		Data	base Re	ef No.			Tested by			_			
					Hosp no	0		_	-		Sam	ple No.			_		Date Tes	ted		-		
Produc				Lot No.		Produc				Lot No.		Produ				Lot		Product		Lot No.		
	Alsevers	ers R144 3466 ed in Alsevers R154 3466				Panel in CellStab Panel Papainised in CellStab				R143 34			in CellN	Media sed in Ce	llModia	1	3 3466 3 3466	Panel in			R146 3466	
ancir	Rh	M	N	S	s	P1	Lua	K	k	Kp ^a	Lea	Leb	Fya	Fyb	Jka	Jk ^b	Other		T EIAT	T	Т	
1	R ₁ wR ₁	0	+	0	+	3	0	0	+	0	0	+	0	+	0	+	Other	3	0		+	+-+
2	R ₁ R ₁	+	0	+	0	0	0	+	+	0	0	0	+	0	+	0		3	0		+	
3	R ₂ R ₂	+	0	0	+	2	0	0	+	0	+	0	+	0	0	+	Cob+, Bgb-	_	3	_		
4	r'r	0	+	0	+	2	0	0	+	0	+	0	+	+	+	0		3	2	A		
5	r"r	0	+	+	0	1	0	0	+	0	0	+	0	+	+	0		3	3			
6	rr	+	+	0	+	1	0	+	0	0	0	+	0	+	+	0		3	3			
7	rr	0	+	0	+	1	0	+	+	0	+	0	0	+	0	+						
8	rr	0	+	0	+	2	0	0	+	+	0	+	+	0	+	0						
9	rr	+	0	+	0	0	0	0	+	0	+	0	+	0	0	+	Cob+					
10	rr	+	0	0	+	4	+	0	+	0	0	+	0	+	0	+						
Auto										_								0				
				-		-															-	\vdash
		+	_		-					-									-		-	+-+
		-				-				-				-			7		-		-	-
		-	-	-	-	-				-				-				_	-		+	+
A matile a	d. Titus		L						L				L	Group				Pheno				
Dilution:	dy Titre	Т			Т		П					Cell	14	Gioup	,			riieiid	,			
Anti-		+	-+	_	+	+		\rightarrow	+	+	+-	Cell	iu	Conc	usion							
Archive		+	_		+			-	+	_	+-	+				Ned	ative r	eacti	ons	in e	nzv	me
Anti		+	_		+		\vdash	\dashv	\neg	+	1	1				_	positiv				_	
Archive					\top				\top		1	1					•					
DAT B	atch No						<u></u>										sible c		mult	ipie	•	
PS		IgG		IgA		lgM		C3c		C30	t I	Ct	1		í	anti	bodies	?				
					$\neg \uparrow$									Repor	ting BM	S:			Checkir	ng BM	S	
Reagen										T			Г		T							
Batch																						

G:\009 DD&R\003 Red Cell Immunohaematology\001 Everyone\Reagent product profiles


No's

Pipette batch No's



Interesting Case 1- Results

- This lady has Anti-Inb, Anti-c and previously Anti-E
- Samples sent to the International Blood Group Reference Laboratory (IBGRL) for exclusions.
- What other tests would you do or consider?
- Further samples?
- Risk of HDFN?
- Blood supply and advice?





Interesting Case 2

- 40 year old pregnant woman
- Last seen by RCI in 2012
- Grouped as O RhD positive by IH1000
- EDD: May 2021 (as well!)



Interesting Case 2- Serology Blood and Transplant

Pipette batch No's


Sample	Name	DA-	TIE	NT	2	Reques	tor		_	_		Data	base R	ef No.		_		Tested by		_		
ate of	Birth	, ,	_			Hosp ne	0			_		Sam	ple No.			_		Date Teste	d	_		
Product Lot No.						Product Lot No.							Product				No.	Product		Lot No.		
anel ir	Alsevers			R144 34		Panel in	CellSta			R143 3	466		in CellA	Media		R16	3 3466	Panel in L			R146 3466	
anel P	apainised	ed in Alsevers R154 3466			Panel Papainised in CellStab				R153 3	466	Panel	Panel Papanised in C			R17	3 3466	EXPIRY D	ATE: 2020.	.11.19			
	Rh	м	N	s	s	P1	Luª	K	k	Kpª	Lea	Le ^b	Fyª	Fy ^b	Jka	Jkb	Other	BRIAT	EIAT	LAT		
1	R ₁ ^w R ₁	0	+	0	+	3	0	0	+	0	0	+	0	+	0	+		4	4	4		
2	R ₁ R ₁	+	0	+	0	0	0	+	+	0	0	0	+	0	+	0		4	3	4		
3	R ₂ R ₂	+	0	0	+	2	0	0	+	0	+	0	+	0	0	+	Cob+, Bgb-	4	4	4		
4	r'r	0	+	0	+	2	0	0	+	0	+	0	+	+	+	0		4	4	4		
5	r"r	0	+	+	0	1	0	0	+	0	0	+	0	+	+	0		4	4	4		
6	rr	+	+	0	+	1	0	+	0	0	0	+	0	+	+	0		4	4	4		
7	rr	0	+	0	+	1	0	+	+	0	+	0	0	+	0	+		4		4		
8	rr	0	+	0	+	2	0	0	+	+	0	+	+	0	+	0		4		4		
9	rr	+	0	+	0	0	0	0	+	0	+	0	+	0	0	+	Cob+	4		4		
10	rr	+	0	0	+	4	+	0	+	0	0	+	0	+	0	+		4		4		
Auto																		0		0		
ntibo	dy Titre		Accession		-			W						Group	o			Pheno				
lution:							T		T			Cell	ld									
nti-														Conc	lusion							
rchive																						
nti									T													
rchive																						
AT B	atch No			***************************************	-				-													
PS		IgG		IgA		IgM		СЗс		C30	d	Ct	ı									
														Repor	ting BMS	S:		C	hecking	BMS		
eagent	:												Π									
atch																						




Interesting Case 2- Results

- Patient is Bombay (Oh)
- Antibody titre of 128
- Risk of HDFN?
- Consultant asks us to do an IgG and IgM titre of the antibody due to the high titre- Why?

IgG ~150 kDa







To summarise...

- Be specific in test requesting!
- Call us to let us know when sending urgent samples or FMH's
- Preparation is key when it comes to patients with known rare antibodies in order for us to be able to meet the supply of blood!
- We are here to offer advice whether it's transfusion advice or sending of samples!



Thank you for listening!

